

# CIRCUMFIX<sup>®</sup> Circumcision Dressing Description

CIRCUMFIX<sup>®</sup> is a ready-made sterile dressing for male circumcision. It is composed of 3 parts. The first part which contacted to the wound is a stretch retention dressing coated with low-adherent substance to minimize pain during undressing. The middle part is high stretch retention dressing which provides suitable compression and absorb exudate. The last part, adhesive tape, helps to fix itself for convenient use.



# Intended use

CIRCUMFIX<sup>®</sup> is used for reduce wound adherence, absorb exudate and provide suitable retention for circumcision wound. CIRCUMFIX<sup>®</sup> can be used by both healthcare and non-healthcare professionals.

# Warning

- Do not use CIRCUMFIX<sup>®</sup>, if the packaging is damaged and/or open prior to use.
- 2. Single-use only.

# Precaution

 Sterility is guaranteed, unless pouch is damaged and/or open prior to use. Do not use CIRCUMFIX<sup>®</sup> if the packaging is damaged.



#### 2. Do not re-sterile

- Do not over compress CIRCUMFIX<sup>®</sup> to the wound, avoiding blood circulated obstruction.
- Should monitor the wound regularly. If there are any swelling, redness, warmth or increased in pain, stop using immediately.

# Contraindication

Do not use the CIRCUMFIX<sup>®</sup>, if the patient has any allergic history to petrolatum, fatty alcohol or any components of the product.

# Instruction for use

- 1. Clean the wound area as per standard protocol.
- 2. Open pouch and pull CIRCUMFIX® out of the tray.
- 3. Remove the cover from the low-adherent part.
- 4. Gently, cover the wound with low-adherent part, then stretch and wrap the remainder of dressing around the wound area, with appropriate tension, until meet the adhesive tape.
- Remove the backing of adhesive tape, then stick it with the wrapped dressing.
- CIRCUMFIX<sup>®</sup> should be changed daily in order to monitor wound characteristics.

# Available sizes

Dressing size	Ref No.	Packing
Infant/Child: 1.5 x 25 cm	CF01	20 pcs
Adult: 2.5 x 40 cm	CF02	20 pcs

# Storage

Keep at room temperature.



# Expiry date

2 years after manufacturing date.

# Symbols used

Symbols	Definitions
8	Do not use if the package is damaged
2	Do not reuse
STERILE R	Sterilized using irradiation
LOT	Batch code or lot number
M	Date of manufacture
	Use by or expiry date
REF	Catalogue number or reference number
	Consult instruction for use
Ť	Keep dry
	Caution, see IFU important cautionary information such as warning and precaution
25 °C	Upper limit temperature at 25°C
e e e e e e e e e e e e e e e e e e e	Manufacturer/ Name and address
EC REP	Authorized representative
CE	Certification mark, CE

## Manufactured by

Novamedic Co., Ltd.

125-125/1 Moo 2 Phraeksamai, Muang, Samutprakarn, 10280, Thailand Tel +662 138 0478

## **Distributed by**

Novatec Healthcare Co., Ltd.

599/15 Moo 5 Srinakarin road, Samrong-Nua, Muang, Samutprakan, 10270 Thailand Tel +662 349 4100 www.novatec.co.th.

## **EU Representative**

mdi Europa GmbH Langenhagener Str. 71 D-30855 Langenhagen Germany Tel. +49-511-3908-95311 www.mdi-europa.com